Fox-Flow Forte 10×

koncentrovaný roztok

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| | | |
|-----|--------------|-------|
| REF | 100-05000-03 | 0,5 L |
| REF | 100-05000-04 | 1 L |

1. Účel použití

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REF 100-05000-04 1 L

1. Intended Use

Fox-Flow Forte 10× concentrated solution is intended for preparation of the sheath fluid for use in flow cytometers. The product is intended for professional laboratory use only.

2. Method Principle

To accurately measure optical properties of cells with a flow cytometer, cells have to pass through the uniformly bright centre of focused laser beams. Light collection optics is focused on the intersection point of cells with the laser beams to pick up fluorescence and scattered light from cells. One purpose of the fluidic system is to move the cells one by one precisely through the sensing region in a liquid stream in such a way that each cell is illuminated by the same amount of light from the lasers [1]. For precise positioning of cells in a liquid jet, the hydrodynamic focusing technique is used in most cytometers [2].

The product Fox-Flow Forte 10× is developed for the preparation of the sheath fluid which complies with the principle described above. The sheath fluid prepared in such a way exhibits low level of autofluorescence, thereby ensuring a high fluorescent signal to noise ratio.

3. Main Components

Sodium chloride, potassium chloride, disodium phosphate, monopotassium phosphate, 2-phenoxyethan-1-ol, sodium fluoride, water, ProClin™ 300. The product does not contain sodium azide.

6. References

- Cossarizza A. et al. Guidelines for the Use of Flow Cytometry and Cell Sorting in Immunological Studies (Second Edition). *Eur J Immunol* 49:1457–1973 (2019).
- Crosland-Taylor P. A Device for Counting Small Particles suspended in a Fluid through a Tube. *Nature* 171: 37–38 (1953).

7. Symbols

REF Catalogue Number

LOT Batch Number

Expiry Date

Temperature Limit

This Way Up

Fragile

IVD In Vitro Diagnostic Medical Device

CE Marking of Conformity

Consult Instruction for Use

Manufacturer

UDI Unique Device Identifier

8. Manufacturer

I.T.A.-Intertact s.r.o.
Pařížská 67/11
110 00 Praha 1
Czech Republic
www.ita-intertact.com
ita@ita-intertact.com

9. Revision History

-

The product is available in two variants:

| Name | Volume | Catalogue No. |
|---------------------------|--------|---------------|
| Fox-Flow Forte 10× (0,5L) | 0,5 L | 100-05000-03 |
| Fox-Flow Forte 10× (1L) | 1 L | 100-05000-04 |

The Safety Data Sheet is available to download on website of the manufacturer.

4. Storage and Handling

Store at 2–25 °C. Use the solution at temperatures specified in an instrument manual. Expiry date of the unopened product is 36 months. After opening, the solution is stable for 60 days at 2–25°C. Do not freeze.

The batch number and expiry date are indicated on the packaging.

5. Procedure

- To obtain 5 L / 10 L of the sheath fluid, dilute Fox-Flow Forte 10× with distilled/ultrapure water in a 1 : 9 (v/v) ratio.
- Mix thoroughly the concentrate Fox-Flow Forte 10×.
- I. To obtain 5 L of the sheath fluid, transfer 4,5 L of distilled/ultrapure water to 0,5 L of the concentrate or fill up to a final volume of 5 L.
II. To obtain 10 L of the sheath fluid, transfer 9 L of distilled/ultrapure water to 1 L of the concentrate or fill up to a final volume of 10 L.
- Mix the solution well.

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Version: 1
Date of Issue: 24.06.2024
Revision Date: -
www.ita-intertact.com